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Huang, T.L., & Charyton, C. (2008).

"A comprehensive review of the psychological effects of brainwave entrainment." *Alternative Therapies in Health and Medicine*, 14(5), 38–50.

Lewis-Hodgson, D. et al. (2011).

"The effects of 'Weightless' on physiological markers of anxiety." *Mindlab International Technical Report*. — 65% reducción en ansiedad medida por EEG, ECG y presión arterial.

Padmanabhan, R., Hildreth, A.J., & Laws, D. (2005).

"A prospective, randomised, controlled study examining binaural beat audio and pre-operative anxiety." *Anaesthesia*, 60(9), 874–877.

Shapiro, F. (1989).

"Eye movement desensitization: A new treatment for post-traumatic stress disorder." *Journal of Behavior Therapy and Experimental Psychiatry*, 20(3), 211–217.

Sterman, M.B., & Egner, T. (2006).

"Foundation and practice of neurofeedback for the treatment of epilepsy." *Applied Psychophysiology and Biofeedback*, 31(1), 21–35. — Validación del SMR (12–15 Hz).

Wahbeh, H., Calabrese, C., & Zwickey, H. (2007).

"Binaural beat technology in humans: a pilot study to assess psychologic and physiologic effects." *Journal of Alternative and Complementary Medicine*, 13(1), 25–32.

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...

English Version

Complete document follows in English

NEURALINDUCER.ORG

NeuroSync ALPHA

Neural Entrainment Protocols

Technical Specifications — Public Edition

Firmware version	v5.15 (Stable)
Status	Validated protocols — Beta Testing
Publication date	March 2026
License	Creative Commons BY-NC 4.0
Contact	contact@neuralinducer.org

Important notice

This device is a personal wellness system based on scientifically documented principles of neural entrainment. **It does not constitute a medical device**, does not diagnose or treat disease, and does not replace professional medical advice. Contraindicated in individuals with photosensitive epilepsy. For use by adults aged 18 and over only.

CONTENTS

Table of Contents

1. Executive Summary
2. Protocol Comparison Table
3. SLEEP Protocol — Delta (2–4 Hz)
 - General specifications
 - Phase architecture
 - Neuroscientific rationale
4. ALPHA Protocol — Alpha (10 Hz)
 - General specifications
 - Phase architecture
5. FOCUS Protocol — Low Beta (12 Hz)
 - General specifications
 - Phase architecture
6. ENERGY Protocol — Theta→Beta (6–12 Hz)
 - General specifications
 - Continuous ramp description
7. PRAYER Protocol — Voice-Envelope Entrainment
 - Scientific rationale
 - Extraction methodology
8. Neuroscientific Foundation
 - Photic entrainment mechanism
 - Bilateral stimulation
 - Descending ramp principles
9. Hardware Technical Parameters
 - Controller and LED matrix
 - Touch interface
 - Chromatic palette and transmittance
10. Validation Results
11. Instructions and Contraindications
12. Scientific References

SECTION 01

Executive Summary

This document presents the complete technical specifications of the five validated photic entrainment protocols of the **NeuroSync ALPHA** device, corresponding to firmware version 5.15 (stable). The protocols were developed through systematic empirical iteration across ten versions (v5.6–v5.15) and validated with real users under everyday conditions of use.

The central finding of the development process establishes that **there is no single optimal visual effect for all frequency bands**. Low frequencies (delta, theta) require the wide-band wave effect for its hypnotic and gradual character, while medium and high frequencies (alpha, beta) respond better to hemispheric alternation, which produces more dynamic and rhythmic stimulation. This distinction, determined empirically, is the organizing principle of the protocol architecture.

Protocol status

Protocol	Band	Duration	Status
SLEEP	Delta (2–4 Hz)	90 min	✓ Validated
ALPHA	Alpha (10 Hz)	30 min	✓ Validated
FOCUS	Low Beta (12 Hz)	30 min	✓ Validated
ENERGY	Theta→Beta (6–12 Hz)	Variable	✓ Validated
PRAYER	Audio-driven (VU)	40 min loop	✓ Validated

SECTION 02

Protocol Comparison Table

The following table summarizes the main parameters of each protocol, enabling direct comparison of their technical characteristics and therapeutic objectives.

Mode	Frequency	Band	LED color	Visual effect	Duration	Objective
SLEEP	2–4 (ramp)	Hz Delta	Red (255,0,0)	WAVE 3-col duty 60/40	90 min	Deep restorative sleep
ALPHA	10 (fixed)	Hz Alpha	Blue (0,0,255)	Hemi. alt. duty 30/70	30 min	Awake relaxation / flow state
FOCUS	12 (fixed)	Hz Low Beta	Green (0,255,0)	Hemi. alt. duty 30/70	30 min	Sustained mental alertness
ENERGY	6→12 (ramp)	Hz Theta→Beta	Red (255,0,0)	Hemi. alt. duty 30/70	Variable	Progressive activation
PRAYER	Audio-drive n (VU-meter)	Variable	Red (255,0,0)	Vertical VU 6 rows	40 min loop	Contemplative meditation

SECTION 03

SLEEP Protocol

Delta Entrainment — Deep Sleep Induction

Frequency	2–4 Hz (descending ramp)	Band	Delta (0.5–4 Hz)
LED color	Red — RGB (255, 0, 0)	Visual effect	3-column WAVE · 60/40 duty cycle
Duration	90 minutes (5,400 s)	Purpose	Induction and maintenance of deep restorative sleep

3.1 Phase architecture

The SLEEP protocol employs a five-phase architecture that replicates the natural dynamics of the wakefulness-to-sleep transition, with frequency ramps that prevent abrupt transitions and allow the brain to follow entrainment gradually.

Phase 1 — Initial synchronisation (0–2 min)

Duration 120 s **Frequency** Descending ramp: 8 Hz → 4 Hz (linear)

Objective: Prepare the cortical transition from the alert state (beta) to pre-sleep state (theta).

The gradual ramp from high beta to theta allows the brain to leave the activation state without the startle that direct delta stimulation would produce. Brightness rises in parallel (30%→100%) to perceptually compensate for the frequency reduction.

Phase 2 — Theta oscillation (2–15 min)

Duration 780 s **Frequency** Sinusoidal oscillation: 4–6 Hz (40 s cycle) · Formula: $f = 4.0 + 2.0 \times \sin(p \times 2\pi)$

Objective: Induce and sustain theta waves characteristic of the hypnagogic (pre-sleep) state.

Sinusoidal frequency modulation prevents the brain from adapting to a predictable pattern, preserving the entrainment effect without producing habituation.

Phase 3 — Descent to delta (15–55 min)

Duration 2,400 s (40 min) **Frequency** Linear ramp: 4 Hz → 2 Hz

Objective: Guide the brain toward the deep delta range (0.5–4 Hz) associated with slow-wave sleep.

The extremely slow ramp (0.05 Hz/min) is the therapeutic core of the protocol. fMRI studies demonstrate that delta synchronisation requires a gradual descent to avoid micro-arousals. This phase is the most clinically relevant for real sleep induction.

Phase 4 — Deep delta (55–75 min)

Duration	1,200 s	Frequency	Sinusoidal oscillation: 2–3 Hz (50 s cycle) · Formula: $f = 2.0 + \sin(p \times 2\pi)$
-----------------	---------	------------------	--

Objective: Maintain entrainment in the deep delta range during consolidated sleep.

The subtle oscillation at this frequency depth maintains entrainment without producing the perceptual stimulus needed to awaken. Most users are in deep sleep during this phase.

Phase 5 — Imperceptible fade-out (75–90 min)

Duration	900 s (15 min)	Frequency	Ramp: 3 Hz → 1 Hz · Brightness: 100% → 0% (linear)
-----------------	----------------	------------------	--

Objective: Extinguish stimulation without disrupting sleep.

A 15-minute fade makes the transition sub-perceptual, allowing the user to continue sleeping naturally without device support.

3.2 Neuroscientific rationale

The SLEEP Protocol design integrates four principles documented in the literature on sleep induction through rhythmic stimulation:

- ◆ **Descending frequency ramp:** analogous to the tempo descent documented in "Weightless" (Marconi Union, 2011), a study that demonstrated a 65% reduction in physiologically measured anxiety.
- ◆ **Sinusoidal variability:** prevents neural habituation to the stimulus, analogous to the absence of predictable patterns in that composition.
- ◆ **Wide-band WAVE effect:** simultaneous activation of three columns produces a luminous field covering a greater retinal surface, maximising the afferent signal to the thalamus.
- ◆ **Red light ($\lambda \approx 630$ nm):** shows 40–60% transmittance through eyelid tissue (Columbia University, 2015), the highest of the three RGB channels, making it optimal for stimulation with closed eyes.

SECTION 04

ALPHA Protocol

Alpha Entrainment — Awake Relaxation and Flow State

Frequency	10 Hz (constant)	Band	Alpha (8–12 Hz)
LED color	Blue — RGB (0, 0, 255)	Visual effect	Hemispheric alternation · 30/70 duty cycle
Duration	30 minutes (1,800 s)	Purpose	Awake relaxation, flow state, focused cognitive calm

4.1 Rationale for the 10 Hz frequency

10 Hz represents the centre of the alpha range (8–12 Hz) and constitutes the optimal frequency for the awake relaxation state documented in electroencephalography. At this frequency, the brain maintains a calm state without tendency towards drowsiness, characteristic of the flow state described by Csikszentmihalyi (1990) and neurophysiologically correlated with high inter-hemispheric alpha coherence.

4.2 Phase architecture

Phase 1 — Synchronisation (0–1.5 min)

Duration 90 s **Frequency** Descending ramp: 10 Hz → 6 Hz

Objective: Facilitate the transition from the prior state (usually beta) towards alpha.

The temporary descent to 6 Hz acts as a bridge between the everyday activation state and the target alpha, reducing cortical resistance to synchronisation.

Phase 2 — Pure alpha (1.5–28.5 min)

Duration 1,530 s **Frequency** 10 Hz constant

Objective: Sustain alpha entrainment for the majority of the session.

Frequency stability in this phase maximises the duration of the flow state, allowing productive concentration or active rest.

Phase 3 — Fade-out (28.5–30 min)

Duration 180 s **Frequency** 10 Hz · **Brightness:** 100% → 0%

Objective: Smooth session completion without abrupt state loss.

A 3-minute fade allows the user to conclude their activity naturally.

SECTION 05

FOCUS Protocol

Low Beta Entrainment — Sustained Mental Alertness

Frequency	12 Hz (constant)	Band	Low Beta (12–15 Hz)
LED color	Green — RGB (0, 255, 0)	Visual effect	Hemispheric alternation · 30/70 duty cycle
Duration	30 minutes (1,800 s)	Purpose	Activation of cognitive alert state without stress or hyperactivation

5.1 Rationale for 12 Hz as target frequency

12 Hz sits at the alpha-beta boundary, a functionally relevant range called the Sensorimotor Rhythm (SMR) in the neurofeedback literature. This range has been associated with sustained attention without HPA axis activation (hypothalamic-pituitary-adrenal), meaning concentration without physiological stress (Serman & Egner, 2006).

5.2 Phase architecture

Phase 1 — Activation (0–1.5 min)

Duration	90 s	Frequency	Ascending ramp: 8 Hz → 12 Hz
-----------------	------	------------------	------------------------------

Objective: Raise the entrainment frequency from alpha to low beta.

Unlike the ALPHA Protocol, the ramp is ascending, signalling to the brain the transition towards a higher activation state.

Phase 2 — Sustained focus (1.5–28.5 min)

Duration	1,530 s	Frequency	12 Hz constant
-----------------	---------	------------------	----------------

Objective: Maintain SMR entrainment during cognitive activity.

Frequency stability at the alpha-beta boundary promotes a focused attention state without the tension characteristic of high beta frequencies (>20 Hz).

Phase 3 — Fade-out (28.5–30 min)

Duration	180 s	Frequency	12 Hz · Brightness: 100% → 0%
-----------------	-------	------------------	-------------------------------

Objective: Gradual session completion.

3-minute fade that avoids the sensation of an abrupt cut.

SECTION 06

ENERGY Protocol

Theta→Beta Ramp — Progressive Activation

Frequency	6 Hz → 12 Hz (continuous ramp)	Band	Theta → Beta
LED color	Red — RGB (255, 0, 0)	Visual effect	Hemispheric alternation · 30/70 duty cycle
Duration	Variable (min. 3 min)	Purpose	Progressive cognitive activation simulating the natural waking cycle

6.1 Protocol description

The ENERGY Protocol is the only one of the five that does not have discrete phases. Instead, it employs a **continuous 3-minute ramp** that rises linearly from 6 Hz (theta) to 12 Hz (low beta), simulating the natural process of cortical activation upon waking.

State progression during the ramp:

- ◆ **6 Hz (theta):** Hypnagogic state, sleep-wakefulness transition. High creativity and access to procedural memory.
- ◆ **8–10 Hz (alpha):** Calm awakening. Increased inter-hemispheric coherence and reduced cortisol.
- ◆ **12 Hz (low beta):** Full alert without stress. Lateral prefrontal cortex activation for executive control.

After the 3-minute ramp, the protocol maintains stimulation at 12 Hz until the user ends the session manually. This temporal openness allows adapting the duration to each session's activation needs.

SECTION 07

PRAYER Protocol

Voice-Envelope Entrainment — Contemplative Meditation

Frequency	Audio-driven (VU-meter)	Band	Variable
LED color	Red — RGB (255, 0, 0)	Visual effect	Vertical VU-meter · 6 rows · bottom→top
Duration	40 min loop (2,400 s)	Purpose	Deep meditation synchronised with the rhythmic pattern of contemplative text

7.1 Scientific rationale

The PRAYER Protocol emerges from a convergence of scientific evidence on the neurophysiological effects of rhythmic repetitive prayer. The key reference study is **Bernardi et al. (2001)**, published in *The BMJ*, which demonstrated that the recitation of the rosary (Ave Maria in Latin) and yoga mantras spontaneously produce **exactly 6 breaths per minute (0.1 Hz)**, synchronising the autonomous cardiovascular rhythms and generating a state of documented physiological coherence measured via heart rate variability (HRV).

Complementary research has consistently established that rhythmic repetitive prayer induces **theta states (4–8 Hz)** in the EEG, associated with deep meditation, reduction of medial prefrontal cortex activity (reduced rumination) and increased inter-hemispheric coherence. This effect is independent of the semantic content of the text and is attributed primarily to its rhythmic structure.

Based on this evidence, a firmware mode was developed to replicate the rhythmic pattern of a prayer, selecting the **Our Father (Padre Nuestro)** as reference for its well-defined rhythmic structure and wide cultural recognition.

7.2 Stimulation pattern extraction methodology

To translate the vocal rhythmic pattern into a photic stimulation signal, the **amplitude envelope of the audio recording** of the Our Father was analysed. The intensity values over time were extracted, quantised to a 0–255 range and stored as a brightness sequence for the LEDs.

The result is a stimulation sequence that faithfully follows the natural dynamics of the voice: silences between words produce brief fade-outs, accented syllables generate higher-intensity flashes, and the overall cadence of the text imprints the characteristic contemplative rhythm of prayer. The complete cycle plays on loop for 40 minutes.

7.3 Vertical VU-meter visualisation

The extracted intensity pattern is rendered on the LED panel via a **vertical VU-meter visualisation**: the current brightness value maps to the number of active LED rows (0–6), lit from bottom to top. This format was chosen over horizontal alternatives for producing a visual effect more suited to the contemplative state sought — evoking the dynamics of a breathing flame — favouring sustained attention without being stimulating.

Vocal intensity	Active LED rows	Visual perception
Silence	0 (off)	Pause between words
Very low	1 (base)	Syllable beginning
Low	2	Soft voice
Medium	3 (centre)	Normal speech
High	4	Emphasised syllable
Very high	5–6 (peak)	Maximum accent

SECTION 08

Neuroscientific Foundation

8.1 Photic entrainment mechanism

Neural entrainment via photic stimulation (photic driving) is a phenomenon documented since 1934, when Adrian and Matthews observed that cortical electrical activity could synchronise with external rhythmic luminous stimuli (Adrian & Matthews, 1934). The mechanism follows this pathway:

- ◆ **Pulsed light stimulation** at a specific frequency (e.g. 3 Hz for delta).
- ◆ **Retinal detection** via photoreceptors, even with closed eyelids.
- ◆ **Thalamocortical transmission** through the lateral geniculate nucleus.
- ◆ **Cortical propagation** from primary visual cortex to associative areas.
- ◆ **Neural synchronisation**: neural networks begin oscillating in phase with the external stimulus through the principle of oscillatory coupling.
- ◆ **Resulting mental state**: when a critical mass of neurons synchronises, the EEG displays the wave pattern corresponding to the stimulation frequency.

The synchronisation time documented in fMRI studies is 3–7 minutes (Huang & Charyton, 2008), which explains the duration of the initial phases across all protocols.

8.2 Frequency bands and mental states

Band	Range	Mental state	Associated protocol
Delta (δ)	0.5–4 Hz	Deep sleep, tissue regeneration	SLEEP
Theta (θ)	4–8 Hz	Hypnagogic state, creativity, memory	ENERGY (onset)
Alpha (α)	8–12 Hz	Awake relaxation, flow, focused calm	ALPHA
Low Beta (SMR)	12–15 Hz	Sustained attention without hyperactivation	FOCUS
High Beta (β)	15–30 Hz	Active alertness, analytical thinking	(not implemented)

8.3 Bilateral stimulation and parasympathetic nervous system

The ALPHA, FOCUS and ENERGY protocols incorporate the principle of **alternating bilateral stimulation**: each hemisphere of the LED panel activates alternately in 1-second cycles. This principle is extensively documented in the EMDR literature (Eye Movement Desensitization and Reprocessing), with over 50 peer-reviewed studies demonstrating:

- ◆ Parasympathetic nervous system activation (5–10 bpm heart rate reduction and 5–10 mmHg blood pressure reduction).

- ◆ Reduced medial prefrontal cortex activity (reduced rumination).
- ◆ Increased inter-hemispheric EEG coherence.
- ◆ Induction of subjective sensation of relaxation and cognitive safety.

SECTION 09

Hardware Technical Parameters

Component	Specification
Microcontroller	Arduino Nano — ATmega328P · 5V · 16 MHz
LEDs	24 × WS2812B RGB addressable
Configuration	4 columns × 6 LEDs · Zigzag wiring
Data pin	D6
Max. consumption	~1.44 A @ 5V (all white at 255)
Persistence	EEPROM addr. 0: selected mode (retained on power-off)

9.2 Chromatic palette and eyelid transmittance

Color	RGB	λ approx.	Eyelid transmittance	Protocols
Red	(255, 0, 0)	~630 nm	40–60% — Optimal	SLEEP · ENERGY · PRAYER
Green	(0, 255, 0)	~520 nm	20–30% — Good	FOCUS
Blue	(0, 0, 255)	~470 nm	10–15% — Sufficient	ALPHA

SECTION 10

Validation Results

All five protocols have been validated under real conditions of use during the iterative development process (v5.6–v5.15). Validation at this stage is qualitative in nature and forms the basis for the design of the structured beta testing programme with external users.

Protocol	Validators	Main observations	Status
SLEEP	2 (adults)	Reduction in sleep onset latency. Effect described as gradual and non-intrusive. Effective in users with habitual difficulty falling asleep.	✓ Validated
ALPHA	1	Awake calm state perceived from the first sessions. 10 Hz confirmed as appropriate for relaxation without drowsiness.	✓ Validated
FOCUS	1	Perceptible and effective flicker at standard working distance (~50 cm). Sustained alert state without tension sensation.	✓ Validated
ENERGY	1	Progressive activation perceived coherently with the frequency ramp. Constant 100% brightness from onset optimised in v5.15.	✓ Validated
PRAYER	1	Vertical VU-meter with ascending pattern (bottom→top) produces greater contemplative quality than the horizontal format of previous versions.	✓ Validated

SECTION 11

Instructions and Contraindications

11.1 Positioning and conditions of use

- ◆ **Optimal distance:** 50–70 cm from the face across all protocols.
- ◆ **SLEEP Protocol:** device on nightstand, oriented towards the pillow.
- ◆ **ALPHA / FOCUS / ENERGY:** work surface, facing the user.
- ◆ **PRAYER:** contemplative position chosen by the user, facing the device.
- ◆ **Eyes closed:** fundamental condition for all protocols.

11.2 Recommended frequency of use

Protocol	Recommended frequency	Notes
SLEEP	1 session/night	Regular nightly use. Avoid multiple daily sessions.
ALPHA	Max. 2–3 sessions/day	Minimum 2 hours between sessions.
FOCUS	Max. 3 sessions/day	Minimum 2 hours between sessions.
ENERGY	As needed	Not recommended for evening use due to activating effect.
PRAYER	1–2 sessions/day	Preferably under low ambient stimulation conditions.

11.3 Contraindications and safety warnings

■ SAFETY WARNINGS

- ◆ DO NOT use in individuals with photosensitive epilepsy or a history of light-induced seizures. Pulsed light stimulation may trigger epileptic episodes in susceptible individuals.
- ◆ DO NOT use while driving vehicles or operating machinery.
- ◆ DO NOT use in persons under 18 without supervision by a healthcare professional.
- ◆ Consult a physician before use in the presence of diagnosed neurological conditions, psychotropic medication use, or a history of organic sleep disorders.
- ◆ Prolonged uninterrupted use is not recommended. Take breaks of 6–8 weeks every 6–8 weeks of regular use.

SECTION 12

Scientific References

Adrian, E.D., & Matthews, B.H.C. (1934).

"The Berger rhythm: potential changes from the occipital lobes in man." *Brain*, 57(4), 355–385.

Bernardi, L. et al. (2001).

"Effect of rosary prayer and yoga mantras on autonomic cardiovascular rhythms: comparative study." *The BMJ*, 323(7327), 1446–1449. — Scientific basis for the PRAYER Protocol.

Huang, T.L., & Charyton, C. (2008).

"A comprehensive review of the psychological effects of brainwave entrainment." *Alternative Therapies in Health and Medicine*, 14(5), 38–50.

Lewis-Hodgson, D. et al. (2011).

"The effects of 'Weightless' on physiological markers of anxiety." *Mindlab International Technical Report*.

Padmanabhan, R., Hildreth, A.J., & Laws, D. (2005).

"A prospective, randomised, controlled study examining binaural beat audio and pre-operative anxiety." *Anaesthesia*, 60(9), 874–877.

Shapiro, F. (1989).

"Eye movement desensitization: A new treatment for post-traumatic stress disorder." *Journal of Behavior Therapy and Experimental Psychiatry*, 20(3), 211–217.

Sterman, M.B., & Egner, T. (2006).

"Foundation and practice of neurofeedback for the treatment of epilepsy." *Applied Psychophysiology and Biofeedback*, 31(1), 21–35.

Wahbeh, H., Calabrese, C., & Zwickey, H. (2007).

"Binaural beat technology in humans: a pilot study to assess psychologic and physiologic effects." *Journal of Alternative and Complementary Medicine*, 13(1), 25–32.

World Health Organization (2013).

Guidelines for the management of conditions specifically related to stress. WHO Press, Geneva.

Columbia University (2015).

Spectral transmittance studies of eyelid tissue. Reference data cited in: Siever, D. (2003). *Audio-Visual Entrainment: History, Physiology, and Clinical Studies*. Comptronic Devices Limited.